

## REF: 660 136

### Antipanic magnetic kit for retractable posts

Antipatico magnetic kit for retractable posts. Dimensions: 25 x 75 mm.

#### Instructions:

1. Apply one or two magnet attachments to the tape.
2. Apply the other to the receiving post or to the receiving clip on the wall.
3. Bring the belt closer to the post to join them through magnet.
4. The mechanism guarantees the tensioning of the belt and the same time the immediate release in the event of a collision.



**Designer Club S.R.L. Via 2 Giugno 28, 21022, Azzate (VA) , Italy**

**Designer Club reserves the right to change specifications and dimensions of products without notice, except for errors or omissions. Copyright Designer Club s.r.l. 2023.**